

# INTCO SYNGUARD

## Nitrile Examination Glove (3 mil)

### 9 inches

- Sky Blue color
- No latex protein to cause allergy
- Excellent softness and wearing fitness
- Undifferentated shelf life as normal gloves
- Well suitable for high cleanliness industries



Intco Medical is a high-tech manufacturing company comitted to R&D, production and marketing of medical disposable products and durable equipment with main business covering medical protection, rehabilitation nursing, physical therapy and disposable examination. INTCO Medical has been exported high quality products to more than 120 countries.

Synguard Nitrile Exam Gloves help reduce the risk of spreading infection or disease. They are non-sterile, latex and powder free. This product is an economical and effective choice to protect against unwanted, dangerous and potentially infectious substances. They complete the range of personal protective equipment designed for the vast majority of professional and medical operations.



Size	S	M	L	XL
Item Number	NGPF7001	NGPF7002	NGPF7003	NGPF7004
Lenght (mm)	230	230	230	230
Width (mm)	80±10	95±10	110±10	120±10
Thickness	Cuff : 0.04 mm / Finger : 0.06 mm / Palm : 0.04 mm			
Physical Properties	Tensile strength : minimum 14 MPa			
	Elongation : minimum 500 %			
Pinholes	AQL 1.5			
Lot Split	10 %	40 %	40 %	10 %
Count	100 gloves/box			
Certifications	• Complies with EN455, EN374, ASTM D6319, ASTM F1671, ASTM 5250 and FDA(510K)			

#### Chemo Therapy Permeation

Yes\*(ask sales rep for full report)  
Chemical permeation and chemotherapy drug resistance data will be provided upon request. Users are advised to review the material safety data sheet for the chemicals being used to determine the required level of protection.

#### Storage & Handling

Store in original packaging in a cool, dry and well-ventilated area, away from dust, direct sunlight, moisture, x-ray and excessive heat above 100°F (37°C)

#### Product Features

- Extremely flexible for bare hand feel
- Light weight yet incredibly strong
- Sleek, compact packaging helps us cut space and cost
- Resists a variety of chemicals
- Textured fingers tips for good grip



#### USA

2915 Ogletown Rd Ste 3830, Newark, DE 19713

1-833-460-2027

info@neobexmedical.com  
neobexmedical.com

#### Canada & Global

Licence : 12269  
85 Albert Street,  
Sherbrooke, QC, J1H 1N1