

Medicom SafeBasics

True Fit Thin (2.8 mil)

9 inches

- Nitrile examination gloves
- Sky-Blue
- Textured
- Powder free
- Non sterile



With over 30 years of scientific expertise, manufacturing facilities across North America, Europe and Asia, and a global network of suppliers, Medicom has an expansive portfolio of personal protective equipment (PPE), infection control solutions and patient care products used by millions of professionals around the world every day. With manufacturing operations strategically distributed across three continents — the company has Health Canada and FDA-approved manufacturing facilities in the United States, China, Taiwan, France and now in Canada— Medicom is uniquely positioned to meet the exponential demand for personal protective products like face mask, medical gowns and other personal protective equipment, even in times of crisis such as global pandemics.

The thinnest SafeBasics nitrile glove offers the ideal combination of quality and price (powder free).



Size	S	M	L	XL
Item Number	1185 B	1185 C	1185 D	1185 E
Length (mm)	≥ 220	≥ 230	≥ 230	≥ 230
Width (mm)	70-90	85-105	100-120	110-130
Thickness	2.8 mil			
Physical Properties	Tensile strength : minimum 14 MPa			
	Elongation : minimum 500 %			
Pinholes	AQL 1.5			
Lot Split	10 %	40 %	40 %	10 %
Count	300 gloves/box			
Certifications	• Meets ASTM D6319-19			

Product Features • The unrivaled strength of nitrile delivers enhanced barrier protection and exceptional fit, while providing the tactile sensitivity necessary for the most complex procedures.

Storage & Handling • Store under cool, dry place.
• Avoid excessive heat (Max. 40°C, 104°F)
• Shield open box from direct sunlight and fluorescent light.

Chemo Therapy Permation • No

