

MEDISAFE Medical

Nitrile Examination Glove (4 mil)

9 inches

- Sky Blue
- Ambidextrous
- Latex Free
- Powder Free
- Textured Finger



These gloves are an economical and effective choice against unwanted and hazardous substances, such as potentially infectious liquids and materials. They complete the range of personal protective equipment (PPE) designed for the vast majority of daily medical tasks and professional operations.



Size	S	M	L	XL
Item Number	MGL-N01-S	MGL-N01-M	MGL-N01-L	MGL-N01-XL
Length (mm)	235±5	235±5	235±5	235±5
Width (mm)	86±4	98±4	107±4	115±4
Thickness	4 mil			
Physical Properties	Tensile strength : minimum 14 MPa			
	Elongation : minimum 500 %			
Pinholes	AQL 1.5			
Lot Split	10 %	40 %	40 %	10 %
Count	100 gloves/box			
Certifications	<ul style="list-style-type: none"> • Meets ISO13485 Standard • Meets EN455 and EN374 • Meets FDA 510K and FDA food grade 			

- Product Features**
- Textured fingertips for exceptional grip and sensitivity
 - Not made with natural rubber latex
 - Anti-static properties
 - Low levels of extractables and particles

- Storage & Handling**
- Store in a cool, dry place. Avoid excessive heat (40°C, 104°F).
 - Shield open box from direct sun and fluorescent lighting.

- Chemo Therapy Permeation**
- Yes*
*Chemical permeation and chemotherapy drug resistance data will be provided upon request. Users are advised to review the material safety data.

