

MEDISCA SAFE-SENSE

Nitrile Examination Glove (4 mil)

9 inches

- Sky-Blue
- Improves user and product protection
- Robust design for superior durability
- Polyurethane coating for easy donning and removal
- Beaded cuff to reduce glove roll down and ripping



Medisca is the global leader in pharmaceutical compounding which is the science responsible for the creation of customized medication. Safe-Sense gloves have been tested and authorized to meet the highest standard for the handling and the manipulation of a great extent of pharmaceutical molecules. From the staggering number of product options available to the ever-changing regulatory guidelines, trust Safe-Sense to provide safety and compliance for all medical needs. Safe-Sense is the new gold standard for simple and reliable PPE product solutions, allowing you to perform your work with peace of mind

These high-end gloves provide maximum protection for non-sterile hazardous compounding, administration, disposal and other related chemotherapy activities. They are resistant to puncture, abrasion and snag. These gloves are also very commonly used in biohazardous environments.



Size	S	M	L	XL
Item Number	5742	5743	5744	5745
Length (mm)	240	240	240	240
Width (mm)	80-90	90-100	100-110	110-120
Thickness	Cuff : 0.05mm / Finger : 0.09 mm / Palm : 0.08 mm			
Physical Properties	Tensile strength : minimum 14 MPa			
	Elongation : minimum 500 %			
Pinholes	AQL 1.5			
Lot Split	10 %	40 %	40 %	10 %
Count	100 gloves/box			
Certifications	<ul style="list-style-type: none"> • Meets ASTM D6978, ASTM D5151 and EN455 standards • Manufactured in ISO 13485 certified facility • Meets FDA regulations for food handling applications 			

- Product Features**
- Produced with distilled water for a silky feeling
 - Not made with natural rubber latex
 - Anti-static properties
 - Textured fingertips to provide grip and sensitivity
 - Chlorinated glove for easy donning and removal
 - Beaded cuff to reduce glove roll

- Chemo Therapy Permeation**
- Yes*
Chemical permeation and chemotherapy drug resistance data will be provided upon request. Users are advised to review the material safety data sheet

- Storage & Handling**
- Store in a cool, dry place

