

MEDLINE SensiCare

Nitrile Examination Glove (3.5 mil)

9.5 inches

- Dark Blue
- Textured
- Powder free
- Non Sterile
- Ambidextrous
- With SmartGuard Film Opening



Through a commitment to make healthcare run better with a culture of responsibility and understanding, Medline is the largest privately-held manufacturer and distributor of medical supplies providing products, education, clinical programs, and services across the continuum of care with offices in 20 countries.

Caring for patients while protecting your staff is more critical than ever. Trust Medline SensiCare for affordable high-performance protection for quality patient care as they are brand recognized by all healthcare professionals. Also, these gloves are recognized as soft and comfortable to provide outstanding tactile sensitivity. Furthermore, the SmartGuard barrier film helps to keep gloves in the box, so you take only what you need to reduce waste which is critical to limit shortages.



Size	S	M	L	XL
Item Number	MDS2584	MDS2585	MDS2586	MDS2587
Length (mm)	240	240	240	240
Thickness	3.5 mil			
Physical Properties	Tensile strength : minimum			
	Elongation : minimum			
Pinholes	AQL 1.5			
Lot Split	10 %	40 %	40 %	10 %
Count	250 gloves/box			
Certifications	• Complies with ASTM D6978-05 and USP 800			

Chemo Therapy Permeation

Yes*(ask sales rep for full report)
Chemical permeation and chemotherapy drug resistance data will be provided upon request. Users are advised to review the material safety data sheet for the chemicals being used to determine the required level of protection.

Storage & Handling

Store in original packaging in a cool, dry and well-ventilated area, away from dust, direct sunlight, moisture, x-ray and excessive heat above 100°F (37°C)

Product Features

- Dark blue powder-free gloves have a soft and comfortable feel while providing outstanding tactile sensitivity for touch perception
- Textured fingertips for enhanced grip
- Color-coded sizing helps you quickly identify the right size every time



USA

2915 Ogletown Rd Ste 3830, Newark, DE 19713

1-833-460-2027

info@neobexmedical.com
neobexmedical.com

Canada & Global

Licence : 12269
85 Albert Street, Sherbrooke, QC, J1H 1N1