

# RONCO NE1

Nitrile Examination Glove (3 mil)

9 inches

- Sky Blue color
- Great chemical resistance
- Better flexibility than vinyl gloves
- 100% latex-free and powder-free
- Medical examination grade glove



Ronco NE1 gloves are a range of specially designed comfort fit gloves for optimum protection and comfort. The NE1 gloves have much better elasticity than Vinyl gloves and offer better puncture resistance than latex gloves. These gloves are latex free and powder free making them an ideal solution for

These gloves are ideal for applications that require dexterity. NE1 gloves are medical grade gloves and can be used in pharmacies, dental offices, and veterinary offices and are a good choice for first responders.



focus on quality



Size	S	M	L	XL
Item Number	735S	735M	735L	735XL
Lenght (mm)	230	230	230	230
Width (mm)	85±4	95±4	106±4	115±4
Thickness	3 mil			
Physical Properties	Tensile strength : minimum 18 MPa			
	Elongation : minimum 500 %			
Pinholes	AQL 1.5			
Lot Split	10 %	40 %	40 %	10 %
Count	100 gloves/box			
Certifications	• Meets ASTM D6319 standard			

**Storage & Handling** • Store in a cool, dry place

**Chemo Therapy Permeation** • No

- Product Features**
- Extremely flexible for bare hand feel
  - Light weight yet incredibly strong
  - Sleek, compact packaging helps us cut space and cost
  - Resists a variety of chemicals
  - Textured fingertips for good grip

