

# NITRILE EXAMINATION GLOVES

## Intco Synguard Gloves

Compliant with ASTM D6319 & ASTM D6978

100 gloves / box  
10 boxes / case

Intco Medical is a high-tech manufacturing company committed to R&D, production and marketing of medical disposable products and durable equipment with main business covering medical protection, rehabilitation nursing, physical therapy and disposable examination. INTOCO Medical has been exported high quality products to more than 120 countries.

**Synguard Nitrile Exam Gloves help reduce the risk of spreading infection or disease. They are non-sterile, latex and powder free. This product is an economical and effective choice to protect against unwanted, dangerous and potentially infectious substances. They complete the range of personal protective equipment designed for the vast majority of professional and medical operations. ideal for healthcare workers, first responders, dental offices, beauty treatments and nursing homes.**



Characteristic		Acceptable limit	
Dimension and tolerance	Width (mm) & Weight (g)	S	70 – 90 mm 2.8g ± 0.3
		M	85 – 105 mm 3.0g ± 0.3
		L	100 – 120 mm 3.3g ± 0.3
		XL	110 – 130 mm 3.6g ± 0.3
	Length (mm) / (Inch)	≥ 230 mm / 9.05"	
	Thickness (mm) / (mil)	Finger	0.06 mm / 2.36 mil
		Palm	0.05 mm / 1.97 mil
		Cuff	0.04 mm / 1.57 mil
Physical requirements	Before Aging	Tensile Strength (MPa)	≥ 14
		Ultimate Elongation (%)	≥ 500
	After Aging	Tensile Strength (MPa)	≥ 14
		Ultimate Elongation (%)	≥ 400
Powder residue		<2mg/glove	

<b>Material</b>	Nitrile
<b>Water Test AQL</b>	1.5 AQL
<b>Description</b>	Non-sterile, latex-free, powder-free, ambidextrous. Chlorinated glove for easy donning and removal, beaded cuff to reduce glove roll down and ripping, textured fingertips to provide grip and sensitivity
<b>Sizes</b>	S, M, L, XL
<b>Color</b>	Sky blue
<b>Conservation</b>	Store in a clean, cool and dry place. Avoid excessive heat (Max 38°C).
<b>Expiration date</b>	3 years after manufacturing date
<b>Manufactured in</b>	China
<b>Certifications</b>	Compliant with ASTM D6319, ASTM D6978, FDA 510K, Health Canada & manufactured in ISO2859
<b>Chemotherapy</b>	Yes, meets ASTM D6978 for resistance to permeation by chemotherapy drugs

**NEOBEX**  
**Manufacturer & Wholesaler**



### USA

2915 Ogletown Rd Ste 3830,  
Newark, DE 19713

1-833-460-2027  
Info@neobexmedical.com  
Neobexmedical.com

### Canada & Global

License : 12269  
85 Albert Street, Sherbrooke, QC,  
Canada J1H 1N1

Copyright © Neobex Medical Inc.  
All Rights Reserved