

NITRILE EXAMINATION GLOVES

Aurelia Amazing Gloves

300 gloves / box
10 boxes / case

Compliant with ASTM D6319, ASTM D6978, FDA 510(K), Health Canada

Aurelia® Gloves has been a top provider of examination gloves in the North American market since 2000. Their innovative and trusted Nitrile and Latex gloves are used across the world, from healthcare to industrial. They are proud to say that they are one of very few glove brands that source directly from their own manufacturing facilities, so they can ensure only the best quality and consistent supply/inventory.

These high-end thin purple blue nitrile examination gloves are specifically designed to offer increased comfort, flexibility & sensitivity. They are resistant to puncture, abrasion and snag. These gloves are also very commonly used in biohazardous environments, ideal for healthcare workers, first responders, dental offices, beauty treatments and nursing homes.



Characteristic		Acceptable limit	
Dimension and tolerance	Width (mm) & Weight (g)	XS S M L XL	70 ± 10 mm 80 ± 10 mm 95 ± 10 mm 110 ± 10 mm 115 ± 10 mm
	Length (mm) / (Inch)		≥ 240 mm / 9.5"
	Thickness (mm) / (mil)	Finger Palm Cuff	0.09 mm / 3.54 mil 0.05 mm / 1.97 mil 0.05 mm / 1.97 mil
Physical requirements	Before Aging	Tensile Strength (MPa)	≥ 14
		Ultimate Elongation (%)	≥ 500
	After Aging	Tensile Strength (MPa)	≥ 14
		Ultimate Elongation (%)	≥ 400
Powder residue		<2mg/glove	

Material	Nitrile
Water Test AQL	≤ 1.5 AQL
Description	Non-sterile, latex & powder-free, ambidextrous. Chlorinated glove for easy donning and removal, beaded cuff to reduce glove roll down and ripping, micro textured to provide grip and sensitivity and produced with pure water for a silky feeling
Sizes	XS, S, M, L, XL
Color	Purple blue
Conservation	Store in a clean, cool and dry place. Avoid excessive heat (Max 38°C).
Expiration date	5 years after manufacturing date
Manufactured in	Malaysia
Certifications	Compliant with ASTM D6319, ASTM D6978, FDA 510(K), FDA regulations for food handling applications, Health Canada & manufactured according to ISO 9001:2015 and ISO 13485:2016 Quality Management Systems.
Chemotherapy	Yes, meets ASTM D6978 for permeation resistance by chemotherapy drugs

Size	Product code
XS	1570-92885
S	1570-92886
M	1570-92887
L	1570-92888
XL	1570-92889



Manufacturer & Wholesaler



USA

2915 Ogletown Rd Ste 3830,
Newark, DE 19713

1-833-460-2027
Info@neobexmedical.com
Neobexmedical.com

Canada & Global

Licence LEIM/MDEL : 12269
85 Albert Street, Sherbrooke, QC,
Canada J1H 1N1

Copyright © Neobex Medical Inc.
All Rights Reserved