

# NITRILE EXAMINATION GLOVES

## Neobex Extended Cuff Gloves

100 gloves / box  
10 boxes / case

Compliant with ASTM D6319 & ASTM D6978, FDA 510(K), Health Canada

In the heart of the Covid-19 pandemic, where markets were struggling to meet global demand faced with shortages and supply difficulties, Neobex developed a quality brand of medical gloves to be part of the growing movement to bring safe and quality medical equipment to its users and facilitate access to health care.

This product is a high-end and effective choice to protect against unwanted, dangerous diseases in environments with the highest risk of contagion. These gloves have been specially developed with a 12-inch sleeve to offer additional protection for the harsh hospital environment where it is essential to protect healthcare workers from infection or disease so that they can continue to play their heroic role in complete safety. Ideal for healthcare workers, first responders, dental offices, beauty treatments and nursing homes.



Characteristic		Acceptable limit												
Dimension and tolerance	Width (mm) & Weight (g)	<table border="1"> <tr><td>S</td><td>80 ± 10 mm</td><td>5.0g ± 0.2</td></tr> <tr><td>M</td><td>95 ± 10 mm</td><td>5.5g ± 0.2</td></tr> <tr><td>L</td><td>110 ± 10 mm</td><td>6.0g ± 0.2</td></tr> <tr><td>XL</td><td>120 ± 10 mm</td><td>6.5g ± 0.2</td></tr> </table>	S	80 ± 10 mm	5.0g ± 0.2	M	95 ± 10 mm	5.5g ± 0.2	L	110 ± 10 mm	6.0g ± 0.2	XL	120 ± 10 mm	6.5g ± 0.2
	S	80 ± 10 mm	5.0g ± 0.2											
	M	95 ± 10 mm	5.5g ± 0.2											
L	110 ± 10 mm	6.0g ± 0.2												
XL	120 ± 10 mm	6.5g ± 0.2												
Length (mm) / (Inch)	≥ 290 mm / 11,5"													
Thickness (mm) / (mil)	<table border="1"> <tr><td>Finger</td><td>0.12 mm / 4,72 mil</td></tr> <tr><td>Palm</td><td>0.09 mm / 3,54 mil</td></tr> <tr><td>Cuff</td><td>0.07 mm / 2,76 mil</td></tr> </table>	Finger	0.12 mm / 4,72 mil	Palm	0.09 mm / 3,54 mil	Cuff	0.07 mm / 2,76 mil							
Finger	0.12 mm / 4,72 mil													
Palm	0.09 mm / 3,54 mil													
Cuff	0.07 mm / 2,76 mil													
Physical requirements	Before Aging	Tensile Strength (MPa)	≥ 14											
		Ultimate Elongation (%)	≥ 500											
	After Aging	Tensile Strength (MPa)	≥ 14											
		Ultimate Elongation (%)	≥ 400											
Powder residue		<2mg/glove												

<b>Material</b>	Nitrile
<b>Water Test AQL</b>	≤ 1.5 AQL
<b>Description</b>	Non-sterile, latex-free, powder-free, ambidextrous. Chlorinated glove for easy donning and removal, beaded cuff to reduce glove roll down and ripping, glove fully textured to provide enhanced grip and sensitivity
<b>Sizes</b>	S, M, L, XL
<b>Color</b>	Sky blue
<b>Conservation</b>	Store in a clean, cool and dry place. Avoid excessive heat (Max 38°C).
<b>Expiration date</b>	3 years after manufacturing date
<b>Manufactured in</b>	China
<b>Certifications</b>	Compliant with ASTM D6319, ASTM D6978, FDA 510(K), Health Canada
<b>Chemotherapy</b>	Yes, meets ASTM D6978 for permeation resistance by chemotherapy drugs

Size	Product code
S	1530-1021000B
M	1530-1021000C
L	1530-1021000D
XL	1530-1021000E



**Manufacturer & Wholesaler**



### USA

2915 Ogletown Rd Ste 3830,  
Newark, DE 19713

1-833-460-2027  
Info@neobexmedical.com  
Neobexmedical.com

### Canada & Global

Licence LEIM/MDEL : 12269  
85 Albert Street, Sherbrooke, QC,  
Canada J1H 1N1

Copyright © Neobex Medical Inc.  
All Rights Reserved