

BLEND VINYL & NITRILE Intco Blue Synmax Exam Gloves

100 gloves / box
10 boxes / case

Compliant with ASTM D5250 FDA 510K & Health Canada

INTCO blend exam gloves of nitrile and vinyl help reduce the risk of spreading infection or disease. They are non-sterile, latex and powder free. This product is an economical and effective choice against unwanted and hazardous substances, such as potentially infectious liquids and materials. They complete the range of personal protective equipment designed for the vast majority of professional and medical operations.



Ideal for healthcare workers, first responders, dental offices, beauty treatments and nursing homes.

Characteristic		Acceptable limit	
Dimension and tolerance	Width (mm) & Weight (g)	S	80 ± 10 mm 4.6g ± 0.2
		M	95 ± 10 mm 5.1g ± 0.2
		L	110 ± 10 mm 5.6g ± 0.2
		XL	120 ± 10 mm 6.1g ± 0.2
	Length (mm) / (Inch)	≥ 230 mm / 9.05"	
	Thickness (mm) / (mil)	Finger	0.05 mm / 1.97 mil
		Palm	0.08 mm / 3,15 mil
		Cuff	0.05 mm / 1.97 mil
Physical requirements	Before Aging	Tensile Strength (MPa)	≥ 15
		Ultimate Elongation (%)	≥ 350
	After Aging	Tensile Strength (MPa)	≥ 15
		Ultimate Elongation (%)	≥ 300
Powder residue		<2mg/glove	

Material	Synthetic Poly-Vinyl Chloride & Nitrile
Water Test AQL	≤ 1.5 AQL
Description	Non-sterile, latex-free, powder-free, ambidextrous. Beaded cuff to reduce glove roll down and ripping
Sizes	S, M, L, XL
Color	Blue
Conservation	Store in a clean, cool and dry place. Avoid excessive heat (Max 38°C).
Expiration date	3 years after manufacturing date
Manufactured in	China
Certifications	Compliant with ASTM D5250 FDA 510K & Health Canada
Chemotherapy	No

Size	Product code
S	1560-BMPF3001
M	1560-BMPF3002
L	1560-BMPF3003
XL	1560-BMPF3004

NEOBEX
Manufacturer & Wholesaler



USA

2915 Ogletown Rd Ste 3830,
Newark, DE 19713

1-833-460-2027
Info@neobexmedical.com
Neobexmedical.com

Canada & Global

Licence LEIM/MDEL : 12269
85 Albert Street, Sherbrooke, QC,
Canada J1H 1N1

Copyright © Neobex Medical Inc.
All Rights Reserved