

# NITRILE EXAMINATION GLOVES

## Medisca Safe-Sense Gloves

100 gloves / box  
10 boxes / case

Compliant with ASTM D6319 & ASTM D6978

Medisca is the global leader in pharmaceutical compounding which is the science responsible for the creation of customized medication. Safe-Sense gloves have been tested and authorized to meet the highest standard for handling and manipulating a great extent of pharmaceutical molecules. From the staggering number of product options available to the ever-changing regulatory guidelines, trust Safe-Sense to provide safety and compliance for all medical needs. Safe-Sense is the new gold standard for simple and reliable PPE product solutions, allowing you to perform your work with peace of mind.



These high-end gloves provide maximum protection for non-sterile hazardous compounding, administration, disposal and other related chemotherapy activities. They are resistant to puncture, abrasion and snag. These gloves are also very commonly used in biohazardous environments, ideal for healthcare workers, first responders, dental offices, beauty treatments and nursing homes.

Characteristic	Acceptable limit													
Dimension and tolerance	Width (mm) & Weight (g)	<table border="1"> <tr> <td>S</td> <td>85 ± 5 mm</td> <td>4.0g ± 0.3</td> </tr> <tr> <td>M</td> <td>95 ± 5 mm</td> <td>4.5g ± 0.3</td> </tr> <tr> <td>L</td> <td>105 ± 5 mm</td> <td>5.0g ± 0.3</td> </tr> <tr> <td>XL</td> <td>115 ± 5 mm</td> <td>5.5g ± 0.3</td> </tr> </table>	S	85 ± 5 mm	4.0g ± 0.3	M	95 ± 5 mm	4.5g ± 0.3	L	105 ± 5 mm	5.0g ± 0.3	XL	115 ± 5 mm	5.5g ± 0.3
	S	85 ± 5 mm	4.0g ± 0.3											
	M	95 ± 5 mm	4.5g ± 0.3											
L	105 ± 5 mm	5.0g ± 0.3												
XL	115 ± 5 mm	5.5g ± 0.3												
Length (mm) / (Inch)	≥ 240 mm / 9.5"													
Thickness (mm) / (mil)	<table border="1"> <tr> <td>Finger</td> <td>0.09 mm / 3,54 mil</td> </tr> <tr> <td>Palm</td> <td>0.08 mm / 3,15 mil</td> </tr> <tr> <td>Cuff</td> <td>0.05 mm / 1.97 mil</td> </tr> </table>	Finger	0.09 mm / 3,54 mil	Palm	0.08 mm / 3,15 mil	Cuff	0.05 mm / 1.97 mil							
Finger	0.09 mm / 3,54 mil													
Palm	0.08 mm / 3,15 mil													
Cuff	0.05 mm / 1.97 mil													
Physical requirements	Before Aging	<table border="1"> <tr> <td>Tensile Strength (MPa)</td> <td>≥ 14</td> </tr> <tr> <td>Ultimate Elongation (%)</td> <td>≥ 500</td> </tr> </table>	Tensile Strength (MPa)	≥ 14	Ultimate Elongation (%)	≥ 500								
	Tensile Strength (MPa)	≥ 14												
	Ultimate Elongation (%)	≥ 500												
	After Aging	<table border="1"> <tr> <td>Tensile Strength (MPa)</td> <td>≥ 14</td> </tr> <tr> <td>Ultimate Elongation (%)</td> <td>≥ 400</td> </tr> </table>	Tensile Strength (MPa)	≥ 14	Ultimate Elongation (%)	≥ 400								
Tensile Strength (MPa)	≥ 14													
Ultimate Elongation (%)	≥ 400													
Powder residue	<2mg/glove													

<b>Material</b>	Nitrile
<b>Water Test AQL</b>	1.5 AQL
<b>Description</b>	Non-sterile, latex-free, powder-free, ambidextrous. Chlorinated glove for easy donning and removal, beaded cuff to reduce glove roll down and ripping, textured fingertips to provide grip and sensitivity and produced with pure water for a silky feeling
<b>Sizes</b>	S, M, L, XL
<b>Color</b>	Sky blue
<b>Conservation</b>	Store in a clean, cool and dry place. Avoid excessive heat (Max 38°C).
<b>Expiration date</b>	3 years after manufacturing date
<b>Manufactured in</b>	China
<b>Certifications</b>	Compliant with ASTM D6319, ASTM D6978, FDA regulations for food handling applications, Health Canada & manufactured in ISO 13485 certified facility
<b>Chemotherapy</b>	Yes, meets ASTM D6978 for permeation resistance by chemotherapy drugs

 **NEOBEX**  
**Manufacturer & Wholesaler**



### USA

2915 Ogletown Rd Ste 3830,  
Newark, DE 19713

1-833-460-2027  
Info@neobexmedical.com  
Neobexmedical.com

### Canada & Global

License : 12269  
85 Albert Street, Sherbrooke, QC,  
Canada J1H 1N1

Copyright © Neobex Medical Inc.  
All Rights Reserved